



Sklep dla diabetyka

<https://en.sugarcubes.pl> +48570037175

sklep@sugarcubes.pl



Enlite Guardian protect skinny patches 10 pcs

catalog number: ENLC10

Kod EAN: 5905090428350

Producer: Sugarcubes

Preparation for shipping: 2-3 days

Price

9,00 EUR

Product description

Protecting patches for sensor Enlite and Guardian Medtronic are a bit of variety in everyday life dominated by readings and measurements. We recommend using patches to secure the sensor when you are more active (sport), especially in children. Beautiful and diverse designs and high-quality workmanship are an additional incentive to use them.

Price for box 10 pc.

Protect patch for Enlite / Guardian Medtronic:

- a wide range of variants and colors
- waterproof and breathable,
- protection and aesthetics,
- additional sensor protection during physical activity

Directions for use:

- make sure your sensor is clean and dry - preferably wipe it with a disinfectant swab
- peel off the patch carefully, starting in one corner
- to avoid touching the glue, start sticking the patch on the skin from the side of the peeled-off side, peeling off the backing a little and sticking the patch on at the same time
- after applying, smooth the patch with your finger to make sure that it sticks well



This product has a decorative function and it's not a medical device.