

Frio DUO cooling case - pink

catalog number: FRIO2RS

Kod EAN: 5039345000302

Producer: FRIO

Preparation for shipping: 1 day

Price

18,00 EUR

Product description

Frio DUO - cooling case for two pens or vials of diabetic medicaments.

Dimensions: 19 cm x 9 cm

The FRIO® case is:

- outer bag made of Cambrelle material,
- well wicking moisture inner pouch cotton with stitching that evenly distributes the non-toxic absorbent keeps insulin in a safe temperature range to activate,
- just immerse in cool water does not require freezing ideal for traveling,
- light and comfortable reusable product note: suitable for storing other temperature sensitive medicines in small packages (when storing other medicines, the manufacturer's instructions on storage conditions should be followed).

Instructions - How to use the case:

- insert the inner insert containing crystals in cold water for 8-12 minutes gently spread the crystals along the visible stitches
- after a moment of soaking the crystals will take the form of a gel with palpable fingers
- when the gel fills the insert chamber, pull the case to drip off the water dry with a paper towel and put in a decorative outer bag
- ready to use note: this is an estimated time we suggest removing the case from the water when the crystals turn into gel filling 3/4 the inner liner and leaving it until the dripping water is absorbed

note: each cold water (tap, river, mineral) can be used to activate the case, do not pour boiling water and immerse in salt water

